


Bibliography

- [1] **Baselga et al.** (2015): Oral presentation. San Antonio Breast Cancer Symposium, Abstract S6-01.
- [2] **Thress et al.** (2014): Poster presentation. ESMO, Abstract #1270P.
- [3] **Schmiegel et al.** (2017): *Mol Oncol.* 11(2):208-219.
- [4] **Saunders et al.** (2016): *Annals of Oncology* 27 (6):149-206.
- [5] **Vidal et al.** (2017): *J Clin Oncol* 35 (suppl 4S; Abstract 607).
- [6] **Jeffers et al.** (2013): *Cancer Res.* 73. SY11-02.
- [7] **Thress et al.** (2015): *Lung Cancer.* 90(3):509-15.
- [8] **Jones et al.** (2016): *J Clin Oncol* 34 (suppl; Abstract 11538).
- [9] **Higgins et al.** (2012): *Clin Cancer Res.* 18(12):3462-9.
- [10] **Schadendorf et al.** (2015): *Eur J Cancer.* 51:S685.

BEAMing Digital PCR Technology



For research use only. Not for use in diagnostic procedures.

© Copyright 2017 - Sysmex Inostics GmbH

Sysmex Inostics GmbH
Falkenried 88, 20251 Hamburg, Germany · Phone +49 40 325907-0 · Fax +49 040 325907-6699 · info@sysmex-inostics.com · www.sysmex-inostics.com

Sysmex Inostics, Inc.
1 Sysmex Way, Mundelein, Illinois 60060, USA · Phone Toll-Free: +1-855-BEAM-DNA (232-6362) · info@sysmex-inostics.com · www.sysmex-inostics.com


Sysmex Corporation
1-5-1 Wakonohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan · Phone +81 78 265-0500 · Fax +81 78 265-0524 · www.sysmex.co.jp

BEAMPERFORMER3.02/2017

www.sysmex-inostics.com

Go Beyond Biopsy with Blood

OncoBEAM™ Technology: BEAMing Digital PCR



Sysmex Inostics' BEAMing Digital PCR technology combines emulsion PCR with magnetic beads and flow cytometry for the highly sensitive detection of mutant tumor DNA molecules.

OncoBEAM™ Advantage: Unparalleled Sensitivity with High Multiplexing Capabilities

Blood-based mutation analysis requires a sensitive technology. Our BEAMing technology delivers increased sensitivity compared to other liquid biopsy methods providing reliable molecular information for timely and non-invasive treatment selection and enhanced prediction of resistance to therapy.




Fig. 1 48% of the analyzed mCRC patients showed ctDNA at <1% fraction

- Data shown at SABCS 2015 highlighted the fact that the OncoBEAM™ PIK3CA blood test was a better predictor of therapy response than the standard of care (SOC) tissue test.¹
- In lung cancer patients, eligible for 3rd generation TKIs, the OncoBEAM™ EGFR blood test accurately identified patients who would benefit from T790M targeted therapy even when tissue as not available for treatment decisions.²
- A meta-analysis of data generated with the OncoBEAM™ CRC assay showed that 48% of the analyzed mCRC patients had a mutant allele fraction (MAF) of between 0.02% – 1.00% (see Fig.1).^{3,4,5}

OncoBEAM™: Analytical Sensitivity

Assay (Mutation Number)		Analytical Sensitivity — Cut-Off (% Mutant Fraction)
EGFR (10)	L858R / Del19	0.03
	T790M / C797S	0.04
KRAS (15) / NRAS (18) / BRAF (2) / HRAS (2)		0.03
AKT1 (1) / ESR1 (12) / PIK3CA (9)		0.02
ALK (12) / ROS1 (3)		0.02 – 0.04*
IDH1 (5) IDH2 (5)		0.02 – 0.04*

* Depending on mutation

OncoBEAM™ Snapshot: Published Performance

Cancer	Marker	Stage	Patient No.	Tissue Analysis	Concordance (%)
NSCLC	EGFR	IV	78	SOC	99 ⁶
	KRAS	IV	78	SOC	92 ⁶
	EGFR		38	SOC	95 ⁷
mCRC	Extended RAS	IV	238	SOC	93.3 ⁸
Breast Cancer	PIK3CA	IV	34	BEAMing	100 ⁹
Melanoma	BRAF V600E	IV	42	Sanger	93 ¹⁰

OncoBEAM™: The Gold Standard for Blood-Based Molecular Cancer Testing of ctDNA



2/3 of the top 15 global pharmaceutical oncology companies are repeat customers



35,000+ samples analyzed across multiple cancer types



Development & Commercialization of first blood-based ALL-RAS IVD KIT for mCRC patients

Chosen technology for OPUS, CRYSTAL, TIGER, CORRECT, BELLE, AURA and a multitude of other landmark trials

