

XN-L CHECK™ ASSAY SHEET

For XN-L series

Material	
XN-L Control Level1	

Lot No.	
52191401	

Exp. Date	
2015/11/17	

Parameters		Target	Limit	Barcode
RBC	10 ⁶ /μL	2.37	0.12	
HGB	g/dL	6.2	0.3	
	mmol/L	3.8	0.2	
HCT	%	17.7	1.8	
MCV	fL	74.7	4.5	
MCH	pg	26.2	2.4	
	fmol	1.624	0.146	
MCHC	g/dL	35.0	5.3	
	mmol/L	21.7	3.3	

Parameters		Target	Limit	Barcode
PLT-I	10 ³ /μL	58	23	
RDW-SD	fL	49.2	7.4	
RDW-CV	%	18.0	2.7	
PDW	fL	9.0	3.6	
PCT	%	0.06	0.03	
MPV	fL	10.5	1.6	
P-LCR	%	25.4	21.6	

There are some parameters which is not displayed on IPU according to the instrument.
Please refer to the package insert for the handling of the XN-L CHECK.
Do not leave control blood in the room temperature over an hour.
Please store it in a refrigerator (2-8°C) immediately after use.

XN-L CHECK™ ASSAY SHEET

For XN-L series

Material	
XN-L Control Level1	

Lot No.	
52191401	

Exp. Date	
2015/11/17	

Parameters	Target	Limit	Barcode
WBC-C 10 ³ /μL	2.41	0.24	
WBC-D 10 ³ /μL	2.41	0.24	
NEUT# 10 ³ /μL	1.10	0.33	
LYMPH# 10 ³ /μL	0.68	0.27	
MONO# 10 ³ /μL	0.19	0.17	
EO# 10 ³ /μL	0.29	0.15	
BASO# 10 ³ /μL	0.15	0.12	

Parameters	Target	Limit	Barcode
IG# 10 ³ /μL	0.30	0.09	
NEUT% %	45.5	9.1	
LYMPH% %	28.4	11.4	
MONO% %	8.0	6.4	
EO% %	12.1	6.1	
BASO% %	6.0	4.7	
IG% %	12.3	3.7	

There are some parameters which is not displayed on IPU according to the instrument.
Please refer to the package insert for the handling of the XN-L CHECK.
Do not leave control blood in the room temperature over an hour.
Please store it in a refrigerator (2-8°C) immediately after use.

XN-L CHECK™ ASSAY SHEET

For XN-L series

Material	
XN-L Control Level1	

Lot No.	
52191401	

Exp. Date	
2015/11/17	

Parameters	Target	Limit	Barcode
RET# 10 ⁵ /μL	0.1225	0.0368	
RET% %	5.17	1.03	
RBC-O 10 ⁶ /μL	2.34	0.23	
HFR %	1.7	11.0	
MFR %	14.2	30.0	
LFR %	84.1	38.0	

Parameters	Target	Limit	Barcode
IRF %	15.9	38.0	
RET-He pg	22.4	2.2	
fmoI	1.390	0.139	
PLT-O 10 ³ /μL	59	41	
IPF# 10 ³ /μL	10.4	9.9	
IPF %	17.6	5.0	

There are some parameters which is not displayed on IPU according to the instrument.
Please refer to the package insert for the handling of the XN-L CHECK.
Do not leave control blood in the room temperature over an hour.
Please store it in a refrigerator (2-8°C) immediately after use.

XN-L CHECK™ ASSAY SHEET

For XN-L series

Material	
XN-L Control Level2	

Lot No.	
52191402	

Exp. Date	
2015/11/17	

Parameters		Target	Limit	Barcode
RBC	10 ⁶ /μL	4.51	0.23	
HGB	g/dL	12.9	0.5	
	mmol/L	8.0	0.3	
HCT	%	36.1	3.6	
MCV	fL	80.0	4.8	
MCH	pg	28.6	2.3	
	fmol	1.775	0.142	
MCHC	g/dL	35.7	5.0	
	mmol/L	22.2	3.1	

Parameters		Target	Limit	Barcode
PLT-I	10 ³ /μL	235	35	
RDW-SD	fL	43.9	6.6	
RDW-CV	%	15.0	2.3	
PDW	fL	9.0	2.7	
PCT	%	0.24	0.07	
MPV	fL	10.1	1.5	
P-LCR	%	21.5	10.8	

There are some parameters which is not displayed on IPU according to the instrument.
Please refer to the package insert for the handling of the XN-L CHECK.
Do not leave control blood in the room temperature over an hour.
Please store it in a refrigerator (2-8°C) immediately after use.

XN-L CHECK™ ASSAY SHEET

For XN-L series

Material	
XN-L Control Level2	

Lot No.	
52191402	

Exp. Date	
2015/11/17	

Parameters	Target	Limit	Barcode
WBC-C 10 ³ /μL	6.93	0.62	
WBC-D 10 ³ /μL	6.91	0.55	
NEUT# 10 ³ /μL	3.36	1.01	
LYMPH# 10 ³ /μL	1.63	0.49	
MONO# 10 ³ /μL	0.59	0.53	
EO# 10 ³ /μL	0.89	0.45	
BASO# 10 ³ /μL	0.44	0.34	

Parameters	Target	Limit	Barcode
IG# 10 ³ /μL	0.92	0.28	
NEUT% %	48.7	7.3	
LYMPH% %	23.6	4.7	
MONO% %	8.5	5.1	
EO% %	12.8	6.4	
BASO% %	6.4	5.0	
IG% %	13.4	4.0	

There are some parameters which is not displayed on IPU according to the instrument.
Please refer to the package insert for the handling of the XN-L CHECK.
Do not leave control blood in the room temperature over an hour.
Please store it in a refrigerator (2-8°C) immediately after use.

XN-L CHECK™ ASSAY SHEET

For XN-L series

Material	
XN-L Control Level2	

Lot No.	
52191402	

Exp. Date	
2015/11/17	

Parameters	Target	Limit	Barcode
RET# 10 ⁶ /μL	0.0951	0.0285	
RET% %	2.11	0.42	
RBC-O 10 ⁶ /μL	4.41	0.44	
HFR %	1.8	11.0	
MFR %	14.4	30.0	
LFR %	83.8	38.0	

Parameters	Target	Limit	Barcode
IRF %	16.2	38.0	
RET-He	pg 23.2	2.3	
	fmol 1.440	0.144	
PLT-O 10 ³ /μL	240	72	
IPF# 10 ³ /μL	43.5	32.6	
IPF %	18.1	5.0	

There are some parameters which is not displayed on IPU according to the instrument.
Please refer to the package insert for the handling of the XN-L CHECK.
Do not leave control blood in the room temperature over an hour.
Please store it in a refrigerator (2-8°C) immediately after use.

XN-L CHECK™ ASSAY SHEET

For XN-L series

Material	
XN-L Control Level3	

Lot No.	
52191403	

Exp. Date	
2015/11/17	

Parameters		Target	Limit	Barcode
RBC	10 ⁶ /μL	5.47	0.27	
HGB	g/dL	17.7	0.7	
	mmol/L	11.0	0.4	
HCT	%	47.8	4.8	
MCV	fL	87.3	5.2	
MCH	pg	32.4	2.6	
	fmol	2.008	0.161	
MCHC	g/dL	37.0	5.2	
	mmol/L	23.0	3.2	

Parameters		Target	Limit	Barcode
PLT-I	10 ³ /μL	575	75	
RDW-SD	fL	42.9	6.4	
RDW-CV	%	13.4	2.0	
PDW	fL	9.4	1.9	
PCT	%	0.57	0.14	
MPV	fL	9.8	1.5	
P-LCR	%	20.6	10.3	

There are some parameters which is not displayed on IPU according to the instrument.
Please refer to the package insert for the handling of the XN-L CHECK.
Do not leave control blood in the room temperature over an hour.
Please store it in a refrigerator (2-8°C) immediately after use.

XN-L CHECK™ ASSAY SHEET

For XN-L series

Material	
XN-L Control Level3	

Lot No.	
52191403	

Exp. Date	
2015/11/17	

Parameters	Target	Limit	Barcode
WBC-C 10 ³ /μL	16.32	1.47	
WBC-D 10 ³ /μL	16.25	1.30	
NEUT# 10 ³ /μL	8.46	2.54	
LYMPH# 10 ³ /μL	3.07	0.92	
MONO# 10 ³ /μL	1.41	1.27	
EO# 10 ³ /μL	2.23	1.12	
BASO# 10 ³ /μL	1.08	0.84	

Parameters	Target	Limit	Barcode
IG# 10 ³ /μL	2.33	0.58	
NEUT% %	52.0	7.8	
LYMPH% %	18.9	3.8	
MONO% %	8.7	4.4	
EO% %	13.7	6.9	
BASO% %	6.7	5.2	
IG% %	14.3	3.6	

There are some parameters which is not displayed on IPU according to the instrument.
Please refer to the package insert for the handling of the XN-L CHECK.
Do not leave control blood in the room temperature over an hour.
Please store it in a refrigerator (2-8°C) immediately after use.

XN-L CHECK™ ASSAY SHEET

For XN-L series

Material	
XN-L Control Level3	

Lot No.	
52191403	

Exp. Date	
2015/11/17	

Parameters	Target	Limit	Barcode
RET# 10 ⁶ /μL	0.0491	0.0147	
RET% %	0.90	0.32	
RBC-O 10 ⁶ /μL	5.33	0.53	
HFR %	1.7	11.0	
MFR %	12.8	30.0	
LFR %	85.5	38.0	

Parameters	Target	Limit	Barcode
IRF %	14.5	38.0	
RET-He pg	24.9	2.5	
famol	1.545	0.155	
PLT-O 10 ³ /μL	592	89	
IPF# 10 ³ /μL	109.3	76.5	
IPF %	18.5	5.0	

There are some parameters which is not displayed on IPU according to the instrument.
Please refer to the package insert for the handling of the XN-L CHECK.
Do not leave control blood in the room temperature over an hour.
Please store it in a refrigerator (2-8°C) immediately after use.