

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 60740804

Exp. date 2016-06-05

L1 : Level 1

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	2.29	5.8	3.6	17.7	77.2	25.3	1.57	32.8	20.3	55
LIMIT	0.11	0.2	0.1	1.1	3.9	1.3	0.08	2.0	1.2	22

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	45.8	16.4	9.0	0.06	9.9	18.2
LIMIT	6.9	2.5	3.6	0.03	1.5	15.5

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	2.95	2.90	1.09	1.14	0.30	0.26	0.16
LIMIT	0.30	0.29	0.33	0.46	0.27	0.13	0.12

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 60740805

Exp. date 2016-06-05

L2 : Level 2

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	4.36	12.3	7.6	36.2	83.1	28.2	1.75	34.0	21.1	231
LIMIT	0.22	0.4	0.2	1.8	4.2	1.4	0.09	2.0	1.3	35

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	43.2	14.7	9.3	0.24	10.2	19.4
LIMIT	6.5	2.2	2.8	0.07	1.5	9.7

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	7.08	6.95	2.80	2.40	0.77	0.70	0.41
LIMIT	0.42	0.42	0.84	0.72	0.69	0.35	0.32

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915