

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 62420804

Exp. date 2016-11-20

L1 : Level 1

	RBC	HGB		HCT	MCV	MCH		MCHC		PLT
	10 ⁶ /μL	g/dL	mmol/L	%	fL	pg	fmol	g/dL	mmol/L	10 ³ /μL
TARGET	2.30	5.9	3.7	17.9	77.9	25.7	1.59	33.0	20.5	57
LIMIT	0.12	0.2	0.1	1.1	3.9	1.3	0.08	2.0	1.2	23

	RDW-SD	RDW-CV	PDW	PCT	MPV	P-LCR
	fL	%	fL	%	fL	%
TARGET	38.4	13.9	8.5	0.05	9.3	14.2
LIMIT	5.8	2.1	3.4	0.03	1.4	12.1

	WBC-D	WBC-C	NEUT#	LYMPH#	MONO#	EO#	BASO#
	10 ³ /μL	10 ³ /μL	10 ³ /μL	10 ³ /μL	10 ³ /μL	10 ³ /μL	10 ³ /μL
TARGET	2.98	2.94	1.06	1.19	0.33	0.25	0.15
LIMIT	0.30	0.29	0.32	0.48	0.30	0.13	0.12

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 62420805

Exp. date 2016-11-20

L2 : Level 2

	RBC	HGB		HCT	MCV	MCH		MCHC		PLT
	10 ⁶ /μL	g/dL	mmol/L	%	fL	pg	fmol	g/dL	mmol/L	10 ³ /μL
TARGET	4.40	12.7	7.9	36.7	83.3	28.9	1.79	34.6	21.5	227
LIMIT	0.22	0.4	0.2	1.8	4.2	1.4	0.09	2.1	1.3	34

	RDW-SD	RDW-CV	PDW	PCT	MPV	P-LCR
	fL	%	fL	%	fL	%
TARGET	41.7	14.2	9.1	0.22	9.8	15.7
LIMIT	6.3	2.1	2.7	0.07	1.5	7.9

	WBC-D	WBC-C	NEUT#	LYMPH#	MONO#	EO#	BASO#
	10 ³ /μL	10 ³ /μL	10 ³ /μL	10 ³ /μL	10 ³ /μL	10 ³ /μL	10 ³ /μL
TARGET	7.03	6.94	2.66	2.56	0.77	0.63	0.41
LIMIT	0.42	0.42	0.80	0.77	0.69	0.32	0.32

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915