

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 62980804

Exp. date 2017-01-15

L1 : Level 1

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	2.31	5.8	3.6	17.5	75.8	25.1	1.56	33.1	20.6	58
LIMIT	0.12	0.2	0.1	1.1	3.8	1.3	0.08	2.0	1.2	23

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	43.6	16.0	8.7	0.06	9.6	15.8
LIMIT	6.5	2.4	3.5	0.03	1.4	13.4

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	3.01	2.96	1.13	1.06	0.36	0.29	0.17
LIMIT	0.30	0.30	0.34	0.42	0.32	0.15	0.13

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 62980805

Exp. date 2017-01-15

L2 : Level 2

	RBC 10 ⁶ /μL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /μL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	4.44	12.2	7.6	35.4	79.7	27.5	1.71	34.5	21.4	232
LIMIT	0.22	0.4	0.2	1.8	4.0	1.4	0.09	2.1	1.3	35

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	42.5	15.2	9.1	0.23	9.9	16.7
LIMIT	6.4	2.3	2.7	0.07	1.5	8.4

	WBC-D 10 ³ /μL	WBC-C 10 ³ /μL	NEUT# 10 ³ /μL	LYMPH# 10 ³ /μL	MONO# 10 ³ /μL	EO# 10 ³ /μL	BASO# 10 ³ /μL
TARGET	6.95	6.82	2.75	2.36	0.75	0.68	0.41
LIMIT	0.42	0.41	0.83	0.71	0.68	0.34	0.32

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915