

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 71000804

Exp. date 2017-07-02

L1 : Level 1

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	2.37	5.8	3.6	18.2	77.0	24.5	1.52	31.9	19.8	58
LIMIT	0.12	0.2	0.1	1.1	3.9	1.2	0.08	1.9	1.2	23

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	46.5	16.8	9.6	0.06	9.8	18.6
LIMIT	7.0	2.5	3.8	0.03	1.5	15.8

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	2.91	2.86	1.06	1.14	0.33	0.23	0.15
LIMIT	0.29	0.29	0.32	0.46	0.30	0.12	0.12

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 71000805

Exp. date 2017-07-02

L2 : Level 2

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	4.50	12.4	7.7	36.7	81.5	27.6	1.71	33.8	21.0	239
LIMIT	0.23	0.4	0.2	1.8	4.1	1.4	0.09	2.0	1.3	36

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	44.7	15.5	9.9	0.24	10.2	19.7
LIMIT	6.7	2.3	3.0	0.07	1.5	9.9

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	6.94	6.81	2.83	2.28	0.74	0.68	0.41
LIMIT	0.42	0.41	0.85	0.68	0.67	0.34	0.32

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 71000806

Exp. date 2017-07-02

L3 : Level 3

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	5.41	16.6	10.3	48.2	89.1	30.7	1.90	34.4	21.4	570
LIMIT	0.27	0.5	0.3	2.4	4.5	1.5	0.10	2.1	1.3	74

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	42.5	13.7	10.4	0.60	10.4	21.8
LIMIT	6.4	2.1	2.1	0.15	1.6	10.9

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	16.41	16.07	7.56	4.35	1.65	1.69	1.16
LIMIT	0.98	0.96	2.27	1.31	1.49	0.85	0.90

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915