

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 73240804

Exp. date 2018-02-11

L1 : Level 1

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	2.33	5.8	3.6	18.1	77.7	24.9	1.54	32.0	19.9	58
LIMIT	0.12	0.2	0.1	1.1	3.9	1.2	0.08	1.9	1.2	23

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	45.9	16.6	9.1	0.06	9.6	16.4
LIMIT	6.9	2.5	3.6	0.03	1.4	13.9

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	2.91	2.82	1.06	1.12	0.32	0.25	0.16
LIMIT	0.29	0.28	0.32	0.45	0.29	0.13	0.12

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 73240805

Exp. date 2018-02-11

L2 : Level 2

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	4.47	12.3	7.6	36.6	81.9	27.5	1.71	33.6	20.9	242
LIMIT	0.22	0.4	0.2	1.8	4.1	1.4	0.09	2.0	1.3	36

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	45.6	15.7	9.2	0.24	9.9	17.2
LIMIT	6.8	2.4	2.8	0.07	1.5	8.6

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	6.96	6.74	2.76	2.33	0.79	0.66	0.42
LIMIT	0.42	0.40	0.83	0.70	0.71	0.33	0.33

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915