

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 83510804

Exp. date 2019-03-10

L1 : Level 1

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	2.29	5.7	3.5	17.4	75.8	24.9	1.54	32.8	20.3	55
LIMIT	0.11	0.2	0.1	1.0	3.8	1.2	0.08	2.0	1.2	22

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	45.0	16.6	8.3	0.05	9.4	14.1
LIMIT	6.8	2.5	3.3	0.03	1.4	12.0

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	3.02	2.92	1.08	1.21	0.31	0.26	0.16
LIMIT	0.30	0.29	0.32	0.48	0.28	0.13	0.12

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 83510805

Exp. date 2019-03-10

L2 : Level 2

	RBC	HGB		HCT	MCV	MCH		MCHC		PLT
	10 ⁶ /uL	g/dL	mmol/L	%	fL	pg	fmol	g/dL	mmol/L	10 ³ /uL
TARGET	4.31	12.0	7.4	35.4	82.1	27.8	1.73	33.9	21.0	223
LIMIT	0.22	0.4	0.2	1.8	4.1	1.4	0.09	2.0	1.3	33

	RDW-SD	RDW-CV	PDW	PCT	MPV	P-LCR
	fL	%	fL	%	fL	%
TARGET	47.2	16.2	8.4	0.22	9.6	14.0
LIMIT	7.1	2.4	2.5	0.07	1.4	7.0

	WBC-D	WBC-C	NEUT#	LYMPH#	MONO#	EO#	BASO#
	10 ³ /uL	10 ³ /uL	10 ³ /uL	10 ³ /uL	10 ³ /uL	10 ³ /uL	10 ³ /uL
TARGET	6.88	6.65	2.69	2.33	0.78	0.67	0.41
LIMIT	0.41	0.40	0.81	0.70	0.70	0.34	0.32

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915