

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 90420804

Exp. date 2019-05-05

L1 : Level 1

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	2.28	5.7	3.5	17.2	75.5	25.0	1.55	33.1	20.6	54
LIMIT	0.11	0.2	0.1	1.0	3.8	1.3	0.08	2.0	1.2	22

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	41.2	15.1	8.1	0.05	9.3	13.5
LIMIT	6.2	2.3	3.2	0.03	1.4	11.5

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	2.77	2.68	1.02	1.08	0.28	0.24	0.15
LIMIT	0.28	0.27	0.31	0.43	0.25	0.12	0.12

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 90420805

Exp. date 2019-05-05

L2 : Level 2

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	4.34	12.0	7.4	35.5	81.8	27.6	1.72	33.8	21.0	224
LIMIT	0.22	0.4	0.2	1.8	4.1	1.4	0.09	2.0	1.3	34

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	42.6	14.8	8.4	0.21	9.6	13.8
LIMIT	6.4	2.2	2.5	0.06	1.4	6.9

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	6.87	6.66	2.75	2.31	0.74	0.66	0.41
LIMIT	0.41	0.40	0.83	0.69	0.67	0.33	0.32

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915