











XN CHECK™ BF ASSAY SHEET

For XN-L series

Material	
BF Control Level1	

Lot No.	
83581301	

Exp. Date	
2019/03/17	


Parameters		Target	Limit	Barcode
WBC-BF	10 ³ /μL	0.077	0.027	
RBC-BF	10 ⁶ /μL	0.025	0.009	
MN#	10 ³ /μL	0.022	0.015	
PMN#	10 ³ /μL	0.055	0.039	
MN%	%	28.7	20.1	
PMN%	%	71.3	49.9	
TC-BF#	10 ³ /μL	0.077	0.027	


Please refer to the package insert for the handling of the XN CHECK BF.
Do not leave control blood in the room temperature over an hour.
Please store it in a refrigerator (2-8°C) immediately after use.









XN CHECK™ BF ASSAY SHEET

For XN-L series

Material	
BF Control Level2	

Lot No.	
83581302	

Exp. Date	
2019/03/17	

Parameters		Target	Limit	Barcode
WBC-BF	10 ³ /μL	0.320	0.080	
RBC-BF	10 ⁶ /μL	0.077	0.019	
MN#	10 ³ /μL	0.090	0.054	
PMN#	10 ³ /μL	0.230	0.138	
MN%	%	28.2	16.9	
PMN%	%	71.8	43.1	
TC-BF#	10 ³ /μL	0.320	0.080	

Please refer to the package insert for the handling of the XN CHECK BF.
Do not leave control blood in the room temperature over an hour.
Please store it in a refrigerator (2-8°C) immediately after use.