

< XN CHECK™ BF ASSAY SHEET >

For XN Series

Lot No. 91051301

Exp. date 2019-07-07

L1 : Level 1

	WBC-BF 10 ³ /μL	RBC-BF 10 ⁶ /μL	MN# 10 ³ /μL	PMN# 10 ³ /μL	MN% %	PMN% %	TC-BF# 10 ³ /μL
TARGET	0.076	0.025	0.020	0.056	25.7	74.3	0.076
LIMIT	0.027	0.009	0.014	0.039	18.0	52.0	0.027

Please refer to the package insert for the handling of the XN CHECK BF.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

< XN CHECK™ BF ASSAY SHEET >

For XN Series

Lot No. **91051302**

Exp. date **2019-07-07**

L2 : Level 2

	WBC-BF 10 ³ /μL	RBC-BF 10 ⁶ /μL	MN# 10 ³ /μL	PMN# 10 ³ /μL	MN% %	PMN% %	TC-BF# 10 ³ /μL
TARGET	0.314	0.072	0.078	0.236	24.9	75.1	0.314
LIMIT	0.079	0.018	0.047	0.142	14.9	45.1	0.079

Please refer to the package insert for the handling of the XN CHECK BF.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.